

Public test session report

BambuLab Standard Indigo Blue PETG CF Blue

MaterialID MAT0170



3DPiceland Labs

Generated with 3DPiceland Engineering Platform v60.0.6

WORKFLOW-RELIABILITY - Bug Fixes and Workflow

Reliability

2026-08-25 19:04:35

Traceability boundary: canonical per-module measured dates are shown below. Session identifier, time of day, operator, printer/slicer profile and environmental conditions are Not recorded and are never inferred.

Measurement provenance

Canonical SQLite measured dates use ISO format; missing dates remain Not recorded.

Module	Measured date
Tensile	Not recorded
Impact	Not recorded
Stiffness	Not recorded

Status
Complete

Result modules
3/3

Specimen/result
records
38

Public details
Approved

Specimen and result quality

Module	Orientation	Average	Std. dev.	CV	Specimens	Confidence	Validation
Tensile	Upright	4,597 MPa	2,757 MPa	0,6%	9	9/10	OK
Tensile	Flat	29,435 MPa	6,919 MPa	0,235%	8	8/10	OK

Module	Orientation	Average	Std. dev.	CV	Specimens	Confidence	Validation
Impact	Upright	89,818 kJ/m ²	29,642 kJ/m ²	0,33%	10	10/10	OK
Impact	Flat	483,779 kJ/m ²	133,906 kJ/m ²	0,277%	10	10/10	OK
Stiffness	Three-point bend	1604,454 MPa	n/a	n/a	1	Complete	OK

Recorded native inputs

Module	Input set	Recorded values
Tensile	Upright force inputs (N)	14, 18, 41, 20, 11, 61, 24, 26, 54,
Tensile	Flat force inputs (N)	199, 198, 202, 222, 214, 202, 82, 212, ,
Impact	Upright needle inputs (%)	5, 4,5, 6,5, 7,5, 4, 3, 7,5, 3,5, 7,5, 4,5
Impact	Flat needle inputs (%)	22, 32, 15, 29, 32, 39, 27, 41, 29, 20
Stiffness	Revolutions / degrees	3 / 70

Approved test notes

Method and equipment context

Purpose-built tensile, pendulum impact and three-point bend apparatus are used with the governed Settings Manager constants consumed by ResultsService. Results are comparative and do not claim accredited-laboratory equipment.

[Testing methodology](#) · [Engineering methodology whitepaper](#) · [Download PDF](#)

reports/test-sessions/MAT0170/